

# -40°C Deep Freezer

Large capacity, excellent deep temperature performance, 8 ergonomic designs

Integrative Design of Evaporator and Shelves Faster Temperature Dropping

This product is designed to store vaccines, blood plasma, biological materials etc. Applications could be found in the electronic and chemical laboratories epidemic prevention, hospitals, animal husbandry areas, blood banks and research institutes

## ● Performance

- Temperature adjustable range:-20°C--40°C
- Microprocessor control, digital display of temperature
- Imported high efficiency Compressor, CFC-free Refrigerant
- Germany EBM brand Fan, high density insulation foam
- Evaporator and Condenser system, guarantee freezing efficiency and uniformity of temperature
- Ergonomic design, more reliable and easy to operate
- High density insulation foam

## ● Safety System

- Multiple alarm functions: Higher or lower temperature alarm, sensor malfunction alarm, power off alarm
- Two kinds of alarm indicators : Audible buzzing alarm, visual flashing alarm
- Start up delay protection

## ● Ergonomic Design



- Double doors design with lock respectively
- Ingenious design of lock, double padlock can be used, providing more safety and reliable
- Door handle, easier to open and close
- Drawer, more convenient to use
- Broad voltage tolerance design, applicable for voltage range 198v-238v
- LED digital display for easy observation
- Equipped with lockable casters, easy to move

# -40°C Deep Freezer



## Specification

| Item Code                         | DW-40L508     |
|-----------------------------------|---------------|
| Capacity(L)                       | 490           |
| Voltage/Frequency                 | 220/50 220/60 |
| Power(W)                          | 430           |
| Temp.range(°C)                    | -20--40       |
| Big/Small Drawer                  | 4/6           |
| Cabinet                           | Upright       |
| Interior Dimensions(mm)(W*D*H)    | 685*650*1228  |
| Exterior Dimensions(mm)(W*D*H)    | 895*840*1860  |
| Packing Dimensions(mm)(W*D*H)     | 965*945*2040  |
| Net/gross Weight(kg)              | 164/170       |
| Container Load(sets) 20'/40'/40'H | 12/24/24      |
| Certificate                       | CE            |



# -40°C Deep Freezer

# -40°C Deep Freezer



| Item Code                         | DW-40L92    | DW-40L188    | DW-40L262     | DW-40L428     | DW-40L626     | DW-40W100   | DW-40W255    | DW-40W380     |
|-----------------------------------|-------------|--------------|---------------|---------------|---------------|-------------|--------------|---------------|
| Capacity(L)                       | 92          | 188          | 262           | 428           | 626           | 100         | 255          | 380           |
| Voltage/Frequency                 | 220/50      | 220/50       | 220/50 220/60 | 220/50        | 220/50        | 220/50      | 220/50       | 220/50 220/60 |
| Power(W)                          | 220         | 290          | 340           | 1200          | 700           | 260         | 410          | 550           |
| Temp.range(°C)                    | -10~-40     | -10~-40      | -10~-40       | -10~-40       | -20~-45       | -20~-40     | -20~-40      | -20~-40       |
| Drawer/Shelve                     | 3/-         | 7/-          | 7/-           | 5/-           | /3            | /           | /            | /             |
| Cabinet                           | Upright     | Upright      | Upright       | Upright       | Upright       | Chest       | Chest        | Chest         |
| Interior Dimensions(mm)(W*D*H)    | 465*440*645 | 380*420*1225 | 480*465*1430  | 760*630*1300  | 760*630*1310  | 500*385*625 | 1040*430*605 | 1376*457*625  |
| Exterior Dimensions(mm)(W*D*H)    | 597*610*860 | 597*610*1475 | 657*685*1665  | 910*835*1990  | 1035*900*1980 | 680*600*805 | 1243*633*838 | 1550*633*838  |
| Packing Dimensions(mm)(W*D*H)     | 659*699*934 | 659*699*1548 | 732*745*1756  | 1065*975*2125 | 1072*954*2130 | 735*645*855 | 1310*750*880 | 1630*755*880  |
| Net/gross Weight(kg)              | 46/51       | 63/71        | 88/93         | 298/314       | 298/314       | 43/46       | 69/81        | 82/89         |
| Container Load(sets) 20'/40'/40'H | 54/108/108  | 27/54/54     | 24/48/48      | 24/48/48      | 12/24/24      | 54/108/162  | 24/54/81     | 18/42/63      |
| Certificate                       | /           | /            | CE            | CE            | /             | /           | CE           | CE            |