



Key Features

- Drawer Structure
- Strong Refrigeration
- Two Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l92.html

Specifications

Technical Data	
Type	-40°C Deep Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	40
Performance	
Cooling performance(°C)	-40
Temperature Range(°C)	-10~-40
Control	
Controller	Microprocessor
Display	LED
Electrical Data	
Power Supply(V/Hz)	220-240/50
Power(W)	220
Electrical Current(A)	2.1
Dimensions	
Capacity(L/Cu.Ft)	92/3.3
Net/Gross Weight(approx)	46/51(kg) 101.4/112.4(lbs)
Interior Dimensions(W*D*H)	465*440*645(mm) 18.3*17.3*25.4(in)

Exterior Dimensions(W*D*H)	597*610*860(mm) 23.5*24.0*33.9(in)
Packing Dimensions(W*D*H)	659*699*934(mm) 25.9*27.5*36.8(in)
Container load (20'/40'/40'H)	54/108/108
Functions	
High/Low Temperature	Y
Sensor Error	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	N/A
Drawers	3
Temperature Recorder	N/A
USB Interface	N/A
CO2 Backup System	
LN2 Backup System	
Chart Recorder	
Others	
Certificate	CE